



Outer capsid protein P8 Antibody

Product Code	CSB-PA801062XA01REN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q85449
Immunogen	Recombinant Rice dwarf virus (isolate Akita) (RDV) Outer capsid protein P8 protein
Raised In	Rabbit
Species Reactivity	Rice dwarf virus (isolate Akita) (RDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice dwarf virus (isolate Akita) (RDV)
Lead Time	Made-to-order (14-16 weeks)