



Outer capsid protein P8 Antibody

Product Code	CSB-PA771187XA01FWE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q85439
Immunogen	Recombinant Rice dwarf virus (isolate Fujian) (RDV) Outer capsid protein P8 protein
Raised In	Rabbit
Species Reactivity	Rice dwarf virus (isolate Fujian) (RDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice dwarf virus (isolate Fujian) (RDV)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.