



Minor outer capsid protein P9 Antibody

Product Code	CSB-PA326854XA01RMJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29078
Immunogen	Recombinant Rice gall dwarf virus (RGDV) Minor outer capsid protein P9 protein
Raised In	Rabbit
Species Reactivity	Rice gall dwarf virus (RGDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice gall dwarf virus (RGDV)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.