



Polyprotein P3 Antibody

Product Code	CSB-PA326724XA01RER
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27502
Immunogen	Recombinant Rice tungro bacilliform virus (isolate Philippines) (RTBV) Polyprotein P3 protein
Raised In	Rabbit
Species Reactivity	Rice tungro bacilliform virus (isolate Philippines) (RTBV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice tungro bacilliform virus (isolate Philippines) (RTBV)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.