



pc1 Antibody

Product Code	CSB-PA769597XA01REQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q85431
Immunogen	Recombinant Rice stripe virus (isolate T) (RSV) pc1 protein
Raised In	Rabbit
Species Reactivity	Rice stripe virus (isolate T) (RSV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice stripe virus (isolate T) (RSV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pc1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.