



S10 Antibody

Product Code	CSB-PA323163XA01REM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19898
Immunogen	Recombinant Rice black streaked dwarf virus (RBSDV) S10 protein
Raised In	Rabbit
Species Reactivity	Rice black streaked dwarf virus (RBSDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice black streaked dwarf virus (RBSDV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	S10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.