



6 Antibody

Product Code	CSB-PA524215XA01RML
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O70791
Immunogen	Recombinant Rice yellow stunt virus (RYSV) (Rice transitory yellowing virus) 6 protein
Raised In	Rabbit
Species Reactivity	Rice yellow stunt virus (RYSV) (Rice transitory yellowing virus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice yellow stunt virus (RYSV) (Rice transitory yellowing virus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.