



p3 Antibody

Product Code	CSB-PA727775XA01FYJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q67897
Immunogen	Recombinant Rice hoja blanca virus (strain cr) (RHBV) (Rice hoja blanca virus (strain Costa Rica)) p3 protein
Raised In	Rabbit
Species Reactivity	Rice hoja blanca virus (strain cr) (RHBV) (Rice hoja blanca virus (strain Costa Rica))
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice hoja blanca virus (strain cr) (RHBV) (Rice hoja blanca virus (strain Costa Rica))
Lead Time	Made-to-order (14-16 weeks)
Target Names	p3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.