



p3 Antibody

Product Code	CSB-PA944894XA01REP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01210
Immunogen	Recombinant Rice stripe virus (isolate M) (RSV) p3 protein
Raised In	Rabbit
Species Reactivity	Rice stripe virus (isolate M) (RSV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rice stripe virus (isolate M) (RSV)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	p3