




gag-pol Antibody

Product Code	CSB-PA361030XA01MHK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03355
Immunogen	Recombinant Moloney murine leukemia virus (isolate Shinnick) (MoMLV) POL protein (560-631aa)
Species Reactivity	Moloney murine leukemia virus (isolate Shinnick) (MoMLV)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Moloney murine leukemia virus (isolate Shinnick) (MoMLV)
Target Names	gag-pol

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: gag-pol antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 23.2 kDa

Observed band size: 26 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.