



bel1 Antibody

Product Code	CSB-PA768245XA01FYM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q87042
Immunogen	Recombinant Simian foamy virus (isolate chimpanzee) (SFVcpz) bel1 protein
Raised In	Rabbit
Species Reactivity	Simian foamy virus (isolate chimpanzee) (SFVcpz)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Simian foamy virus (isolate chimpanzee) (SFVcpz)
Lead Time	Made-to-order (14-16 weeks)
Target Names	bel1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.