



bel1 Antibody

Product Code	CSB-PA526336XA01FEK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O56862
Immunogen	Recombinant Feline foamy virus (FFV) (Feline syncytial virus) bel1 protein
Raised In	Rabbit
Species Reactivity	Feline foamy virus (FFV) (Feline syncytial virus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Feline foamy virus (FFV) (Feline syncytial virus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	bel1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.