



bel1 Antibody

Product Code	CSB-PA320363XA01HQF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14353
Immunogen	Recombinant Human spumaretrovirus (SFVcpz(hu)) (Human foamy virus) bel1 protein
Raised In	Rabbit
Species Reactivity	Human spumaretrovirus (SFVcpz(hu)) (Human foamy virus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Human spumaretrovirus (SFVcpz(hu)) (Human foamy virus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	bel1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.