



env Antibody

Product Code	CSB-PA333080XA01BKG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25506
Immunogen	Recombinant Bovine leukemia virus (isolate Belgium LB285) (BLV) env protein
Raised In	Rabbit
Species Reactivity	Bovine leukemia virus (isolate Belgium LB285) (BLV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bovine leukemia virus (isolate Belgium LB285) (BLV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	env
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.