



env Antibody

Product Code	CSB-PA325800XA01EFY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16082
Immunogen	Recombinant Equine infectious anemia virus (strain WSU5) (EIAV) env protein
Raised In	Rabbit
Species Reactivity	Equine infectious anemia virus (strain WSU5) (EIAV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Equine infectious anemia virus (strain WSU5) (EIAV)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	env