



env Antibody

Product Code	CSB-PA361281XA01KAF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04502
Immunogen	Recombinant Kirsten murine leukemia virus (KiMSV) env protein
Raised In	Rabbit
Species Reactivity	Kirsten murine leukemia virus (KiMSV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Kirsten murine leukemia virus (KiMSV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	env
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.