



env Antibody

Product Code	CSB-PA319877XA01MIT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10259
Immunogen	Recombinant Mouse mammary tumor virus (strain BR6) (MMTV) env protein
Raised In	Rabbit
Species Reactivity	Mouse mammary tumor virus (strain BR6) (MMTV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mouse mammary tumor virus (strain BR6) (MMTV)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	env