



env Antibody

Product Code	CSB-PA311626XA01SHF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02837
Immunogen	Recombinant Simian immunodeficiency virus agm.grivet (isolate AGM gr-1) (SIV-agm.gri) (Simian immunodeficiency virus African green monkey grivet) env protein
Raised In	Rabbit
Species Reactivity	Simian immunodeficiency virus agm.grivet (isolate AGM gr-1) (SIV-agm.gri) (Simian immunodeficiency virus African green monkey grivet)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Simian immunodeficiency virus agm.grivet (isolate AGM gr-1) (SIV-agm.gri) (Simian immunodeficiency virus African green monkey grivet)
Lead Time	Made-to-order (14-16 weeks)
Target Names	env
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.