



env Antibody

Product Code	CSB-PA339356XA01ATC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33498
Immunogen	Recombinant Avian retrovirus RPL30 env protein
Raised In	Rabbit
Species Reactivity	Avian retrovirus RPL30
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian retrovirus RPL30
Lead Time	Made-to-order (14-16 weeks)
Target Names	env
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.