



env Antibody

Product Code	CSB-PA326550XA01BKE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25505
Immunogen	Recombinant Bovine leukemia virus (isolate American VDM) (BLV) env protein
Raised In	Rabbit
Species Reactivity	Bovine leukemia virus (isolate American VDM) (BLV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bovine leukemia virus (isolate American VDM) (BLV)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	env