



env Antibody

Product Code	CSB-PA356035XA01UAM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03398
Immunogen	Recombinant UR2 sarcoma virus (UR2SV) (Avian sarcoma virus (strain UR2)) env protein
Raised In	Rabbit
Species Reactivity	UR2 sarcoma virus (UR2SV) (Avian sarcoma virus (strain UR2))
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	UR2 sarcoma virus (UR2SV) (Avian sarcoma virus (strain UR2))
Lead Time	Made-to-order (14-16 weeks)
Target Names	env
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.