



V-FOS Antibody

Product Code	CSB-PA360489XA01FAC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01102
Immunogen	Recombinant FBJ murine osteosarcoma virus (FBJ-MSV) (Finkel-Biskis-Jinkins murine osteosarcoma virus) V-FOS protein
Raised In	Rabbit
Species Reactivity	FBJ murine osteosarcoma virus (FBJ-MSV) (Finkel-Biskis-Jinkins murine osteosarcoma virus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	FBJ murine osteosarcoma virus (FBJ-MSV) (Finkel-Biskis-Jinkins murine osteosarcoma virus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	V-FOS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.