



gag Antibody

Product Code	CSB-PA324194XA01HHZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21435
Immunogen	Recombinant Hortulanus murine leukemia virus (HoMuLV) (Mus hortulanus virus) gag protein
Raised In	Rabbit
Species Reactivity	Hortulanus murine leukemia virus (HoMuLV) (Mus hortulanus virus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Hortulanus murine leukemia virus (HoMuLV) (Mus hortulanus virus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gag
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.