



# gag Antibody

Product Code	CSB-PA361026XA01ATN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03342
Immunogen	Recombinant Avian spleen necrosis virus gag protein
Raised In	Rabbit
Species Reactivity	Avian spleen necrosis virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian spleen necrosis virus
Lead Time	Made-to-order (14-16 weeks)
Gene Names	gag