



gag Antibody

Product Code	CSB-PA326864XA01MKC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29167
Immunogen	Recombinant Murine leukemia virus (strain BM5 eco) gag protein
Raised In	Rabbit
Species Reactivity	Murine leukemia virus (strain BM5 eco)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Murine leukemia virus (strain BM5 eco)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gag
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.