



gag Antibody

Product Code	CSB-PA320848XA01RAM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11269
Immunogen	Recombinant Radiation murine leukemia virus gag protein
Raised In	Rabbit
Species Reactivity	Radiation murine leukemia virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Radiation murine leukemia virus
Lead Time	Made-to-order (14-16 weeks)
Gene Names	gag