



gag Antibody

Product Code	CSB-PA365910XA01MHL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03334
Immunogen	Recombinant Moloney murine sarcoma virus (MoMSV) gag protein
Raised In	Rabbit
Species Reactivity	Moloney murine sarcoma virus (MoMSV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Moloney murine sarcoma virus (MoMSV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gag
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.