



V-FPS Antibody

Product Code	CSB-PA360551XA01ATI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00541
Immunogen	Recombinant Avian sarcoma virus (strain PRCII) V-FPS protein
Raised In	Rabbit
Species Reactivity	Avian sarcoma virus (strain PRCII)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian sarcoma virus (strain PRCII)
Lead Time	Made-to-order (14-16 weeks)
Target Names	V-FPS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.