



gag Antibody

Product Code	CSB-PA365922XA01ARF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03373
Immunogen	Recombinant Avian erythroblastosis virus (strain ES4) gag protein
Raised In	Rabbit
Species Reactivity	Avian erythroblastosis virus (strain ES4)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian erythroblastosis virus (strain ES4)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gag
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.