



gag Antibody

Product Code	CSB-PA301701XA01EFZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69732
Immunogen	Recombinant Equine infectious anemia virus (strain Wyoming) (EIAV) gag protein
Raised In	Rabbit
Species Reactivity	Equine infectious anemia virus (strain Wyoming) (EIAV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Equine infectious anemia virus (strain Wyoming) (EIAV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gag
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.