



gag Antibody

Product Code	CSB-PA333196XA01FDB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26806
Immunogen	Recombinant Friend murine leukemia virus (isolate FB29) (FrMLV) gag protein
Raised In	Rabbit
Species Reactivity	Friend murine leukemia virus (isolate FB29) (FrMLV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Friend murine leukemia virus (isolate FB29) (FrMLV)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	gag