



gag Antibody

Product Code	CSB-PA341338XA01MHO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32594
Immunogen	Recombinant Moloney murine sarcoma virus (strain ts110) (MoMSV) gag protein
Raised In	Rabbit
Species Reactivity	Moloney murine sarcoma virus (strain ts110) (MoMSV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Moloney murine sarcoma virus (strain ts110) (MoMSV)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	gag