



gag Antibody

Product Code	CSB-PA365909XA01ATM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03327
Immunogen	Recombinant Avian sarcoma virus (strain Y73) gag protein
Raised In	Rabbit
Species Reactivity	Avian sarcoma virus (strain Y73)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian sarcoma virus (strain Y73)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	gag