




HBZ Antibody

Product Code	CSB-PA314490XA01HQG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C745
Immunogen	Recombinant Human T-cell leukemia virus 1 (isolate Caribbea HS-35 subtype A) (HTLV-1) HBZ protein (1-209aa)
Species Reactivity	Human T-cell leukemia virus 1 (isolate Caribbea HS-35 subtype A) (HTLV-1)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Human T-cell leukemia virus 1 (isolate Caribbea HS-35 subtype A) (HTLV-1)
Target Names	HBZ

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: HBZ Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 29 kDa

Observed band size: 36 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.