



JUN Antibody

Product Code	CSB-PA361514XA01ATE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05411
Immunogen	Recombinant Avian sarcoma virus (strain 17) JUN protein
Raised In	Rabbit
Species Reactivity	Avian sarcoma virus (strain 17)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian sarcoma virus (strain 17)
Lead Time	Made-to-order (14-16 weeks)
Target Names	JUN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.