



MYC Antibody

Product Code	CSB-PA319189XA01ASR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10395
Immunogen	Recombinant Avian myelocytomatosis virus CMII MYC protein
Raised In	Rabbit
Species Reactivity	Avian myelocytomatosis virus CMII
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian myelocytomatosis virus CMII
Lead Time	Made-to-order (14-16 weeks)
Target Names	MYC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.