



MYC Antibody

Product Code	CSB-PA356768XA01ASS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06295
Immunogen	Recombinant Avian myelocytomatosis virus HBI MYC protein
Raised In	Rabbit
Species Reactivity	Avian myelocytomatosis virus HBI
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian myelocytomatosis virus HBI
Lead Time	Made-to-order (14-16 weeks)
Gene Names	MYC