



# MYC Antibody

Product Code	CSB-PA319420XA01ATB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12523
Immunogen	Recombinant Avian retrovirus OK10 MYC protein
Raised In	Rabbit
Species Reactivity	Avian retrovirus OK10
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian retrovirus OK10
Lead Time	Made-to-order (14-16 weeks)
Gene Names	MYC