



MYC Antibody

Product Code	CSB-PA325084XA01FAT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21438
Immunogen	Recombinant Feline leukemia virus FTT MYC protein
Raised In	Rabbit
Species Reactivity	Feline leukemia virus FTT
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Feline leukemia virus FTT
Lead Time	Made-to-order (14-16 weeks)
Target Names	MYC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.