



MYC Antibody

Product Code	CSB-PA365521XA01AST
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01110
Immunogen	Recombinant Avian myelocytomatosis virus MC29 MYC protein
Raised In	Rabbit
Species Reactivity	Avian myelocytomatosis virus MC29
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian myelocytomatosis virus MC29
Lead Time	Made-to-order (14-16 weeks)
Gene Names	MYC