



orfC Antibody

Product Code	CSB-PA803424XA01WBH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q88936
Immunogen	Recombinant Walleye dermal sarcoma virus (WDSV) orfC protein
Raised In	Rabbit
Species Reactivity	Walleye dermal sarcoma virus (WDSV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Walleye dermal sarcoma virus (WDSV)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	orfC