



pol Antibody

Product Code	CSB-PA361031XA01ASW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03360
Immunogen	Recombinant Avian reticuloendotheliosis virus pol protein
Raised In	Rabbit
Species Reactivity	Avian reticuloendotheliosis virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Avian reticuloendotheliosis virus
Lead Time	Made-to-order (14-16 weeks)
Target Names	pol
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.