



pol Antibody

Product Code	CSB-PA320733XA01BAB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10272
Immunogen	Recombinant Baboon endogenous virus (strain M7) pol protein
Raised In	Rabbit
Species Reactivity	Baboon endogenous virus (strain M7)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Baboon endogenous virus (strain M7)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	pol