



gag-pol Antibody

Product Code	CSB-PA323124XA01BKB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19560
Immunogen	Recombinant Bovine immunodeficiency virus (strain R29) (BIV) (Bovine immunodeficiency-like virus) gag-pol protein
Raised In	Rabbit
Species Reactivity	Bovine immunodeficiency virus (strain R29) (BIV) (Bovine immunodeficiency-like virus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bovine immunodeficiency virus (strain R29) (BIV) (Bovine immunodeficiency-like virus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gag-pol
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.