



pol Antibody

Product Code	CSB-PA329106XA01FAH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31792
Immunogen	Recombinant Feline endogenous virus ECE1 pol protein
Raised In	Rabbit
Species Reactivity	Feline endogenous virus ECE1
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Feline endogenous virus ECE1
Lead Time	Made-to-order (14-16 weeks)
Gene Names	pol