



pol Antibody

Product Code	CSB-PA323900XA01FAM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19028
Immunogen	Recombinant Feline immunodeficiency virus (strain San Diego) (FIV) pol protein
Raised In	Rabbit
Species Reactivity	Feline immunodeficiency virus (strain San Diego) (FIV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Feline immunodeficiency virus (strain San Diego) (FIV)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	pol