



pol Antibody

Product Code	CSB-PA330116XA01RAN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31795
Immunogen	Recombinant Radiation murine leukemia virus (strain Kaplan) pol protein
Raised In	Rabbit
Species Reactivity	Radiation murine leukemia virus (strain Kaplan)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Radiation murine leukemia virus (strain Kaplan)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pol
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.