



pol Antibody

Product Code	CSB-PA328654XA01FDB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26809
Immunogen	Recombinant Friend murine leukemia virus (isolate FB29) (FrMLV) pol protein
Raised In	Rabbit
Species Reactivity	Friend murine leukemia virus (isolate FB29) (FrMLV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Friend murine leukemia virus (isolate FB29) (FrMLV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pol
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.