



pol Antibody

Product Code	CSB-PA336111XA01MAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35956
Immunogen	Recombinant Maedi visna virus (strain KV1772) (MVV) (Visna lentivirus) pol protein
Raised In	Rabbit
Species Reactivity	Maedi visna virus (strain KV1772) (MVV) (Visna lentivirus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Maedi visna virus (strain KV1772) (MVV) (Visna lentivirus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pol
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.