



pol Antibody

Product Code	CSB-PA356023XA01AGG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03356
Immunogen	Recombinant AKV murine leukemia virus (AKR (endogenous) murine leukemia virus) pol protein
Raised In	Rabbit
Species Reactivity	AKV murine leukemia virus (AKR (endogenous) murine leukemia virus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	AKV murine leukemia virus (AKR (endogenous) murine leukemia virus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pol
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.